


Shelf Liners



LINER DIMENSIONS	A	B

- Tough high density polyethylene protects cabinets and shelving from acid spills
- Over a dozen sizes available to fit most standard size shelves
- Custom sizes available
- Liners are 1/8" thick.
- Polypropylene available on special order

Standard Dimensions

Length (A)	Width (B)
18	12
24	12
24	18
30	12
30	18
36	12
36	18
42	12
42	18
42	24
48	12
48	18
48	24
60	12
60	18
60	24

All sizes in inches.

Specifications:

Shelf liners shall conform to ASTM # 1412 for corrosion resistance and shall be manufactured by Scientific Plastics, Inc. Of Kansas City, Ks.

Liners for doors also available.