


ADA Trays



Item Number	A	B	C	D
ADA 2212-201	12"	22"	2"	2"
ADA 2316-201	16"	23"	2"	2"
ADA 2619-201	19"	26"	2"	2"

- Protects cabinets and shelving from acid spills
- Three sizes of trays available
- 1" or 2" depths
- Standard color is natural, special colors available
- 3/16" tray walls for rigidity
- One-piece, rotationally molded polyolefin
- Rounded corners for easy cleaning
- Dishwasher safe

Specifications: Laboratory trays shall be a minimum of 3/16" thick, one piece, stress free, molded of virgin polyolefin resins with smooth interior walls and coved corners for easy cleaning. Laboratory trays and shelf liners shall conform to ASTM #1412 for corrosion resistance and material content and shall be manufactured by Scientific Plastics, Inc. of Kansas City, KS.