




Door Liners

Scientific Plastics door and shelf liners are made from .125 HDPE and manufactured to your specifications. We can cut the liner to accommodate hinges, latches, vent holes, and holes for attaching to the door. Submit drawing for pricing.



Shelf Liners

Scientific can cut a liner to fit your metal or wood shelf with attachment holes. Submit drawing for pricing.



Encapsulated Shelf

Designed to fit our acid cabinet liners, Scientific's encapsulated shelf is manufactured using a steel frame rail and then encapsulating with polyethylene. This eliminates exposed metal or wood shelves which can corrode over time in an acid storage cabinet.