



Counter Top Sink With Splash

One piece molded polyolefin counter top sink features molded in tailpiece, marine edge, splash guard, and 4" deck for water fixture mounting. Commonly used on free standing sink basin Scientific Plastics part number W31085-000. Standard color is black. Other colors are available.



Part Number	A	B	C	D	E	F	G	H	I	J	K	OUTLET
W31082-000	17-3/8"	26-3/4"	9-3/4"	12-3/4"	26-1/2"	17-1/2"	4"	5"	3"	8"	32-5/8"	1-1/2"